




**DEUTSCH REMOVABLE BRIDGE ZR® APPOINTMENT PROTOCOL (6 implant protocol)**

**Appointment #1**

**Implant level impression with implant open tray copings  
Opposing impression**

**Appointment #2**

**Verified Impression - Placement of verification jig, luting of jig, and impression below jig,  
and over the jig with provided custom tray**



### Appointment #3


Centric and Vertical Relationship Records (bite registration) and photos (at eye level with patient, straight head on shot at rest, smile, and high smile positions).

### Appointment #4

Wax Try-in with screw retained baseplate, check protrusive and excursive balance, phonetics, and esthetic. Photos


### Appointment #5

DRB ZR® Abutment Try-in and torque to 30 NcM.



XRAY verify complete abutment seating.

Place secondary gold telescopic copings using numbers facing the buccal.



Ensure that the Tertiary structure seats passively over the gold secondary telescopes (if binding occurs, adjust tertiary until passive), and that the tissue stops seat on tissue.



Ensure that the gold copings are NOT rebounding on gingival mucosa, and that they are completely seated and DRY.

Using the provided cement, place enough cement in all tertiary receptacles in one pass, and covering the bottom of the receptacle surface completely with cement.

Place the tertiary over the secondary gold telescopes and seat completely (should be just shy of the gold margins). Excess cement will vent through the buccal vents.

Ensure that no cement vents around the abutment itself, and have an assistant use a micro brush to clean any excess in the event that cement leaks around the abutment margin while holding the tertiary structure in place.

The cement can be light cured to tac in place through the vents, and DO NOT remove until the final cement cure time of **8 minutes**.

Remove the tertiary structure, and **DO NOT** attempt to reseat on the working cast.



**Remove and place DRB ZR® abutments back on cast in the correct orientation.**

**Send all back to Deutsch Dental Arts, LLC.**

**Appointment #6 (OPTIONAL)**

**Full hardware try in with final wax try in**

**Appointment #7**

**Final delivery and seating of prosthesis.**

**Place and torque all DRB ZR® abutments to 30 NcM**

**Carefully place DRB prosthesis and seat completely.**

**Ensure proper protrusive and excursive balance and phonetics.**

**Cover DRB ZR® abutment screw access holes with teflon tape or silicone, and seal with flowable composite ensuring a concave shape which does NOT rise above the access hole, and does NOT obstruct proper seating of the gold secondary telescopes in any way.**

**Instruct patient on removal and observe the patient perform this function at least 3 times, and correct hygiene protocols prior to patient dismissal.**

Deutsch Dental Arts, LLC  
12425 W. Bell Road Suite 103  
Surprise, AZ 85378  
623.972.2648 *office*  
623.236.5249 *cell*  
Email : [arian@deutschlab.com](mailto:arian@deutschlab.com)

